

Neuro-Sponge (Cottonoid) Patties — Fabco-Equivalent

Latex-free rayon neurosurgical patties with locator string, sold sterile, 10 patties per card / 1 card per pouch / 20 pouches per box. Direct-from-manufacturer supply via our Authorized-Distributor relationship with American Surgical Company (Salem, MA).

RFQ # _____ Date _____

Bill To

Ship To

Ray-Cot®

Rayon Neurosurgical Pattie

- 12% greater absorbency
- 2× stronger than the leading competitor
- Single barium line
- X-ray detectable locator string

Product Family	Ray-Cot® / Delicot®
Manufacturer	American Surgical Company
Origin	Salem, Massachusetts, USA
Distributor	Sterile Mate LLC (Authorized)
Material	Rayon (latex-free)
Sterility	Sterile
Radiopaque	Single barium sulfate line
Locator String	Yes (X-ray detectable)
Packaging	10/card · 1 card/pouch · 20 pouches/box
Box Count	200 patties / box
Lead Time	Off-the-shelf where available

Clinical Applications

- Neurosurgery — cortical protection, hemostasis, and tissue retraction
- Micro-surgery / ENT / ophthalmology — delicate soft-tissue absorption
- Spine — dural protection and CSF containment
- General surgery — pediatric, oncology, and vascular procedures

Fabco → American Surgical Ray-Cot® Size Cross-Reference

Direct 1:1 replacement across the fourteen most-ordered Fabco pattie sizes. All American Surgical codes below ship 200 patties per box.

Size (in)	Size (mm)	Fabco (typical)	American Surgical Ray-Cot®
-----------	-----------	-----------------	----------------------------

¼ × ¼	6 × 6	✓	60-02
¼ × 1½	6 × 38	✓	60-16
¼ × 3	6 × 76	✓	60-04
½ × ½	13 × 13	✓	60-06
½ × 1	13 × 25	✓	60-07
½ × 1½	13 × 38	✓	60-08
½ × 2	13 × 51	✓	60-29
½ × 3	13 × 76	✓	60-09
¾ × ¾	19 × 19	✓	60-10
1 × 1	25 × 25	✓	60-12
1 × 3	25 × 76	✓	60-13
1½ × 1½	38 × 38	✓	60-14
1½ × 3	38 × 76	✓	60-22
3 × 3	76 × 76	✓	60-15

Also available: Ray-Cot® strips without locator string (Codes 61-01 through 61-08, sizes ¼×6 to 3½×6 in), and the ultra-thin **Delicot®** micro-surgical pattie line for ENT, ophthalmology, and cranial base work.

Quote & Order Instructions

Field	Value
Contract	Non-GPO — competitive vs. Fabco / Codman / Aspen
Unit of Sale	Box (200 patties)
Unit Price	Request current pricing — quoted per size / total-box volume
Payment Terms	Net 30 upon credit approval · Card / ACH / wire accepted
Freight	FOB Salem, MA — freight-included offers available on 5+ box orders
Quote Validity	30 days from download / issue date
Lead Time	Off-the-shelf when possible; allocated orders confirmed at PO

To place this order: Email orders@sterilemate.com or reply to your Sterile Mate account manager with (1) hospital name / bill-to / ship-to, (2) the American Surgical codes and box quantities you need, and (3) your PO number. MBE / HUB / SDVOSB supplier-diversity credit applies where eligible.