



## Exam Indigo Nitrile Disposable Gloves

### Product Code: AINPF

AMMEX<sup>®</sup> Professional Exam Indigo Nitrile Gloves deliver consistent, dependable barrier protection against bodily fluids, bloodborne pathogens, and contaminants. They feature excellent fit, feel, and comfort for easy manipulation and a dexterous grip. The fit and feel rivals latex, without triggering latex sensitivities.

### KEY USER BENEFITS:

- Improved grip and dexterity
- High level of comfort
- Medical grade - meets AQL 1.5
- Color coding with indigo, blue, and black

### FEATURES:

- Textured
- Powder Free
- Latex Free
- Available Sizes: S-XL



Medical



Dental



Tattoo



Child Care



Industrial

### COLOR OPTIONS

3 mil



ABNPF APFN



## GENERAL DESCRIPTION

Color	Indigo	Usage	Single Use
Material	Nitrile	Interior	Powder Free
Grade	Exam	Exterior	Textured
Shape	Ambidextrous	Process	Polymer Coated
Cuff	Beaded	Latex Content	Latex Free
Type	Non-Sterile	Sizes	S-XL

## PHYSICAL PROPERTIES









Size	Unit	Small	Medium	Large	X-Large	Standard Method	
Weight	±0.2 g	-	3.5	-	-	-	
Width	±5 mm	85	95	105	115	ASTM D412 ASTM D3767	
Thickness	Cuff	±1 mil / ±0.02 mm	3 / 0.07	3 / 0.07	3 / 0.07		3 / 0.07
	Palm	±1 mil / ±0.02 mm	3 / 0.08	3 / 0.08	3 / 0.08		3 / 0.08
	Finger	±1 mil / ±0.02 mm	4 / 0.11	4 / 0.11	4 / 0.11	4 / 0.11	
Minimum Length	mm	240					
Maximum Powder Content	mg/glove	2				ASTM D6124	
Maximum Protein Content	µg/dm <sup>2</sup>	N/A				ASTM D5712	

		Before Aging	After Aging	Standard Method
Minimum Tensile Strength	MPa	14	14	ASTM D412
Minimum Ultimate Elongation	%	500	400	

## PRODUCT QUALITY AND REGULATORY COMPLIANCE

QMS Audit Standards	ISO 9001:2015, ISO 13485:2016
Regulatory and Standard Compliance	California Proposition 65; REACH; FDA 21 CFR 177.2600; FDA 510(k); ASTM D6319
AQL Level	1.5

## PRODUCT CODE AND PACKAGING

Size	Small	Medium	Large	X-Large
Product Code	AINPF42100	AINPF44100	AINPF46100	AINPF48100
Box Barcode	 6 97383 95878 5 S Box	 6 97383 95889 1 M Box	 6 97383 95900 3 L Box	 6 97383 95911 9 XL Box
Case Barcode	 6 97383 95944 7 S Case	 6 97383 95955 3 M Case	 6 97383 95966 9 L Case	 6 97383 95977 5 XL Case
Quantity	100 Gloves Per Box, 10 Boxes Per Case			