



## Exam Blue Nitrile Disposable Gloves

### Product Code: GPNHD


Gloveworks® Exam Blue Nitrile Gloves feature an extended cuff—longer than a standard glove—to offer greater protection over your forearm and wrist from bodily fluids, bloodborne pathogens, and contaminants. They are more than twice as thick as standard nitrile gloves for heavy-duty jobs.

#### KEY USER BENEFITS:

- Medical grade - meets AQL 1.5
- Extended cuff provides forearm protection
- 50 count
- Thick for heavy-duty jobs

#### FEATURES

- Textured
- Powder Free
- Latex Free
- Available Sizes: M-XL



### GENERAL DESCRIPTION

Color	Blue	Usage	Single Use
Material	Nitrile	Interior	Powder Free
Grade	Exam	Exterior	Textured
Shape	Ambidextrous	Process	Polymer coated
Cuff	Beaded	Latex Content	Latex Free
Type	Non-Sterile	Sizes	M-XL

### PHYSICAL PROPERTIES







Size	Unit	Medium	Large	X-Large	Standard Method	
Weight	±0.2 g	9.9	-	-	-	
Width	±5 mm	95	105	115	ASTM D412 ASTM D3767	
Thickness	Cuff	±1 mil / ±0.02 mm	6 / 0.14	6 / 0.14		6 / 0.14
	Palm	±1 mil / ±0.02 mm	7 / 0.18	7 / 0.18		7 / 0.18
	Finger	±1 mil / ±0.02 mm	8 / 0.19	8 / 0.19	8 / 0.19	
Minimum Length	mm	295				
Maximum Powder Content	mg/glove	2			ASTM D6124	
Maximum Protein Content	µg/dm <sup>2</sup>	N/A			ASTM D5712	

		Before Aging	After Aging	Standard Method
Minimum Tensile Strength	MPa	14	14	ASTM D412
Minimum Ultimate Elongation	%	500	400	

### PRODUCT QUALITY AND REGULATORY COMPLIANCE

QMS Audit Standards	ISO 9001:2015; ISO 13485:2016
Regulatory and Standard Compliance	California Proposition 65; REACH; FDA 21 CFR 177.2600; FDA 510(k); ASTM D6319
AQL Level	1.5

### PRODUCT CODE AND PACKAGING

Size	Medium	Large	X-Large
Product Code	GPNHD64100	GPNHD66100	GPNHD68100
Box Barcode	 6 97383 40362 9 <b>M Box</b>	 6 97383 40363 6 <b>L Box</b>	 6 97383 40364 3 <b>XL Box</b>
Case Barcode	 6 97383 41362 8 <b>M Case</b>	 6 97383 41363 5 <b>L Case</b>	 6 97383 41364 2 <b>XL Case</b>
Quantity	50 Gloves Per Box, 10 Boxes Per Case		